

# INJECTION *OR* INFECTION

It is important to use disposable needles.  
It is more important to destroy them after use.



**INSTANT**  
**NEEDLE BURNER &**  
**SYRINGE DESTROYER**

Also Available  
**REPLACEABLE  
CUTTER MODEL**

**NULIFE and YOU**  
Together let's save lives



## I N S T A N T NEEDLE BURNER & SYRINGE DESTROYER



A 1999 study by **World Health Organisation** indicates that about 12 billion injections are administered each year. It is shocking but true, infected needles cause an estimated 8-16 million infection of **Hepatitis B**, 2.3 to 4.7 million with **Hepatitis C** and 80,000-1,60,000 with the dreaded **HIV Virus**. All put together, these infections cause an estimated 1.3 million premature deaths.

*Medical professionals are not spared either, every year more than one million medics in the United States alone are affected by these infected needles while handling and recapping causing an unknown number of infections viz: HIV, HBV, HCV, etc.*

**Used, disposable syringes and needles are recycled back into the medical profession by unscrupulous elements, unknown to healthcare professionals/practitioners.**

**Used needles destroyed by mechanical cutters remain infected, therefore they should be destroyed by burning as they pose a threat to even those who are handling them...**

### Then what is the solution ?

It is becoming an accepted practice worldwide to avoid recapping and to destroy the needles and disposable syringes immediately after use on the spot. This is the only way to eliminate any possibility of contamination due to accident or fraud.

Instant Needle Burner & Syringe Destroyer is the safe, permanent and economical way for all concerned healthcare workers to deal with this

### Advantages of **NULIFE DOTS** Instant Needle Burner & Syringe Destroyer:

- ▶ Compact Device that fits easily on the desktop.
- ▶ Portable, safe and cost effective system which creates a safe environment for Doctors, Nurses, Pathologists, Healthcare Professionals and Patients.
- ▶ After use, the contaminated syringe with needle is inserted in the destroyer which generates a high temperature above 1600 degrees centigrade, thereby completely burning the needle and reducing it to ashes in just 2 seconds.
- ▶ Pulling the cutter handle, cuts off the syringe hub and rendering it unuseable, which in turn eliminates chances of it's sales back into unsuspecting markets by unscrupulous elements.
- ▶ Very high grade imported steel cutter is used which does not rust or becomes blunt even after long usage.
- ▶ Removable discharge tray for easy disposal of needle ash and syringe hubs.
- ▶ Shock proof, fuse protection and negligible power consumption as

Size	165mmx115mmx120mm Handy desktop model.
Weight	1.7 kg.
Power	230 voltage±10%, 50 hertz, single phase.
Protection	2 ampere 3/4" cartridge fuse
<b>Also Available Heavy Duty Model</b>	

*We reserve the right to modify/change the specifications/details without notice.*



## **MRK HEALTHCARE**

B4/B5, Byculla Service Industries,  
D. K. Marg, Mumbai 400 027.  
Tel: 91-22-23748371/2/3. Fax: 91-22-23748370.  
E-mail: mrkhealthcare@vsnl.com  
Website: www.mrkhealthcare.com

